

FIRST DRAW & FLUSH LEAD TEST SAMPLING INSTRUCTIONS

STOP PLEASE READ ALL INSTRUCTIONS COMPLETELY BEFORE COLLECTING YOUR WATER SAMPLE STOP

This test kit represents one sampling location. If you require testing on multiple locations you will need to purchase multiple test kits.

1. First Draw sample: Choose a collection time when the water has not been used for at least six (6) hours. Before you use any water from the tap, fill the sample vial labeled "First Draw" as soon as you turn on the tap.
2. Flush sample: Run the water from the tap for 3-5 minutes and then fill the sample vial labeled "Flush".
3. Complete the lab form, ***making sure to indicate the collection date and time you sampled***, type of water sampled (treated/untreated) and sample identification which identifies where you collected your sample (i.e., kitchen sink, well head, etc.). Please be sure to retain a copy of the completed lab form prior to shipment. Return the completed paperwork and water samples to the shipping container.
4. **It is not necessary to ship by an overnight delivery service. Ground, or Saver delivery is sufficient. Cost of shipping samples to the laboratory is at your expense. NTL offers discounted UPS shipping, if interested call for a quote or to obtain your label**

Testing Location: **National Testing Laboratories, Ltd.
556 South Mansfield Street
Ypsilanti, MI 48197**

5. PLEASE NOTE: Testing typically takes approximately 5-7 business days from receipt of samples. Results will be emailed upon completion.

NOTE: The lab reserves the right to reject water samples that appear to be excessively contaminated as they may cause potential damage to laboratory equipment.

If you have any questions regarding this First Draw & Flush Lead Test or any other services we provide, please contact a Technical Service Representative at National Testing Laboratories at 1-800-458-3330 for assistance.

Thank you for choosing National Testing Laboratories, Ltd., for your water testing services!