

## WATERCHECK 1 & 2 SAMPLING INSTRUCTIONS



**PLEASE READ ALL THE INSTRUCTIONS COMPLETELY BEFORE COLLECTING YOUR SAMPLE**



**One test kit represents one sampling location. If you require testing on multiple locations you will need to purchase multiple test kits.**

1. Freeze the blue gel ice pack for at least 8 hours prior to sample collection. Return to cooler just prior to sampling and/or return shipping.
2. Remove any aerators, strainers, mixers or purifiers from the tap.
3. Run the cold water for approximately 3-5 minutes, then adjust flow to a thin stream to avoid splashing or overfilling.
4. Fill the bottle to the neck, cap securely and return to the shipping container.
5. Complete the enclosed lab form, **making sure to indicate the collection date and time you sampled**, type of water (treated/untreated) and sample identification which identifies where you collected your sample (i.e., kitchen sink, well head, etc.). Retain a copy of the lab form for your records and return the completed form to the shipping container.
6. *It is not necessary to ship by an overnight delivery service. Ground, Saver, or 2nd day delivery is sufficient. Cost of shipping samples to the laboratory is at your expense. Samples may be dropped off at the testing location Monday through Friday between 8am - 4pm.*

Testing Location: **National Testing Laboratories, Ltd.  
556 South Mansfield Street  
Ypsilanti, MI 48197**

PLEASE NOTE: Testing typically takes approximately 5-7 business days from receipt of the samples. Reports will be emailed upon completion.

**NOTE: *The lab reserves the right to reject water samples that appear to be excessively contaminated as they may cause potential damage to laboratory equipment.***

If you have any questions regarding this Watercheck 1 & 2 Test or any other services we provide, please contact a Technical Service Representatives at 1-800-458-3330 for assistance

Thank you for choosing National Testing Laboratories, Ltd., for your water testing services!